

Bone Expansion & Collection



BONE EXPANDER KIT SD-BE

Scan the QR code to get more information of the item on this page.

Used for the patient with narrow bone ridge or bad bony tissues.

Unique design of expander drill.



Contouring Bur

·Flattens the slope of the bone.



Saw disk

Thickness
0.25mm



Guide drill



- It is used for single implant.
 - After guide drilling on narrow bone ridge, make sequential expansion using expander drill.
 - It gives initial fixation power.
- Ex) Expanding up to 4.0 for 4.5 size of implant placement.

Expander drill



·Recommended RPM : 40~50

Handpiece condenser



Wrench



Clinical Case

