

# Sinus Lift



## **TOLA KIT** SD-TOLA : Tool of Lateral Approach

Scan the QR code to get more information of the item on this page.

Composition of tools forming sinus lateral window in safe and speedy way.

Used in surgery cases : Residual bone within 1~3mm, Perforated membrane at crestal approach, Placement of multiple implants.



### **Sinus curette**

- Used to lift sinus membrane manually.
- Recommended to use by its sequential order.

#### **\*Detail**




TOLA-01



TOLA-02



TOLA-03



TOLA-04

