

Dental Handpiece



SD-TORQUE SD-TQ

Scan the QR code to get more information of the item on this page.

Used to fix or separate upper structures in implant surgery.

Accessible to the molar area which is hard to implant with wrench.

Installation



· Insert driver into the hole.

Torque control



· Easy torque setting.

Operation



· Turn the handle clockwise or counterclockwise to tighten and loosen.



Head
Contra angle driver head

Latch key

* Autoclave OK

Stable approach and grip

Torque handle

* Drivers are included, (1,2mm)

(Short) 23.5mm

(Long) 30mm

Torque selector

Easy torque setting (5~35N/cm)